



December 28, 2002

Mesia Hovsepian

# 10/072,484

"New Tech Articulator"

Filed 2/7/02

ABSTRACT

The Articulator is used in the dental field. It combines the upper and lower models of the impression taken by the dentist of the patient's bite. When the dental technician receives the impression, he pours the master models until they dry and the models are combined in the articulator with the use of plaster, to register patient's bite to further design the prosthetic.

The New-Tech Articulator is designed to be used without the use of plaster. This will eliminate the use of materials such as, plaster, glue, oil, etc. which are time consuming and messy (as shown in Fig. 4).